

## Delta-Cast® Prints



- Unique printed casting material
- Non-Fibreglass Casting Tape
- Smooth finish and soft edges
- Available in 7 popular patterns

# Delta-Cast® Prints

- Non-Fibreglass Casting Tape
- Easy to mould
- Quick and easy application
- Lightweight cast
- Reliable, consistent working time
- Set time 3 – 5 minutes / full weight-bearing in approximately 20 minutes

## World's only printed casting material

Delta-Cast Prints are non-fibreglass cast tapes which are strong, lightweight and durable. The patterns are colourfast and will not interfere with radiograph viewing. Delta-Cast® Prints are comprised of a patented printed substrate which features extensible yarns that easily shape to body contours, providing enhanced moulding capabilities. Delta-Cast® Prints are non-fibreglass casts that produce less dust during cast removal compared to fibreglass casts.

## Smooth finish and soft edges

The patented non-fibreglass casting material combines strength and rigidity with smooth finish and soft edges for enhanced patient comfort. The smooth surface will not snag clothes and bed linens.

## Available in 7 popular patterns



Paw



Pastel



Teddy  
Bears



Dinosaurs



Football



Camouflage



Christmas\*

## Delta-Cast Prints do not need to be used in conjunction with any other cast tape

Order Number	Pattern	Size	Unit of Measure
72273-08	Football	5cm x 3.6m	1 x 10 rolls
72273-09	Football	7.5cm x 3.6m	1 x 10 rolls
72273-22	Camouflage	5cm x 3.6m	1 x 10 rolls
72273-06	Camouflage	7.5cm x 3.6m	1 x 10 rolls
72273-05	Teddy Bear	5cm x 3.6m	1 x 10 rolls
72273-02	Dinosaur	5cm x 3.6m	1 x 10 rolls
72273-07	Pastel	7.5cm x 3.6m	1 x 10 rolls
72273-20	Paw	5cm x 3.6m	1 x 10 rolls
72273-21	Paw	7.5cm x 3.6m	1 x 10 rolls
72273-03	Christmas Holiday*	5cm x 3.6m	1 x 10 rolls
72273-04	Christmas Holiday*	7.5cm x 3.6m	1 x 10 rolls

\* Only available from September to December.



**LATEX  
FREE**  
FORMULA