

# Delta-Cast® Soft

Semi-Rigid Polyester Cast Tape



The flexible  
solution

# Delta-Cast® Soft

## Semi-Rigid Polyester Cast Tape

Delta-Cast® Soft is a semi-rigid polyester (non-fiberglass) cast tape, which allows for greater range of movement, earlier mobility and reduces some of the problems of rigid immobilization, such as muscle atrophy. It is mainly used in pediatric applications, such as club foot treatment and serial casting, as well as customized support.

### Advanced Technology

Delta-Cast® Soft is a unique, high-tech construction of polyester and polyurethane resin that combines safety, functionality and comfort into one product. Delta-Cast® Soft will hold its form and not delaminate, even as a removable cast.

### Versatility

Delta-Cast® Soft is an excellent choice for functional immobilization, providing support while allowing a limited range of motion. By varying the layers of the cast, the support can be customized to the patient's needs.

### Comfort

Delta-Cast® Soft polyester substrate delivers a smooth surface and soft edges that provide greater patient comfort and satisfaction. A semi-rigid cast minimizes the rehabilitation period by allowing the patient to continue with some daily activities.

### Ease of Application and Removal

The unique Delta-Cast® Soft construction offers consistent unrolling tension, ideal work and set times and a highly conformable, multi-directional knit. Delta-Cast® Soft can be removed with scissors or unwrapping. The use of a cast saw is optional.



Code	Dimensions	Box Contents	Selling Unit
7270802	2.5 cm x 1.8 m	10 rolls	Box
7270812	5 cm x 3.6 m	10 rolls	Box
7270822	7.5 cm x 3.6 m	10 rolls	Box
7270832	10 cm x 3.6 m	10 rolls	Box
7270842	12.5 cm x 3.6 m	10 rolls	Box
7270805	5 cm x 3.6 m	10 rolls	Box
7270815	7.5 cm x 3.6 m	10 rolls	Box
7270825	10 cm x 3.6 m	10 rolls	Box
7270808	5 cm x 3.6 m	10 rolls	Box
7270818	7.5 cm x 3.6 m	10 rolls	Box
7270828	10 cm x 3.6 m	10 rolls	Box
7270801	5 cm x 3.6 m	10 rolls	Box
7270811	7.5 cm x 3.6 m	10 rolls	Box
7270821	10 cm x 3.6 m	10 rolls	Box
7270809	5 cm x 3.6 m	10 rolls	Box
7270819	7.5 cm x 3.6 m	10 rolls	Box



Available in a wide range of attractive colors and unique paw pattern