

Delta-Xpress™ product range

Delta-Xpress™

Delta-Xpress™ offers four different designs to immobilise lower leg, arm, wrist and finger fractures.

Delta-Xpress™ Ordering Information

BSN Code	Size	Splints per pack
Wrist and finger splint		
72713-01	M	1 x 10
72713-02	L	1 x 10

Wrist splint		
72712-01	M	1 x 10
72712-02	L	1 x 10

Lower leg splint		
72711-00	S	1 x 10
72711-01	M	1 x 10
72711-02	L	1 x 10

Elbow splint		
72710-00	S	1 x 10
72710-01	M	1 x 10
72710-02	L	1 x 10

Benefits

Fast, precise and secure immobilisation

- **Fast** - time saving instant application
- **Precise** - ready-to-apply in clinical position
- **Secure** - easy fixation to ensure limb remains in optimal healing position
- **100% radiolucency**
- **Ultra light-weight design**



- **Precise Clinical Positioning**
- **Quick Immobilisation**
- **Cost Effective**

Scan QR Code for more information



Why is it the fastest and most precise **immobilisation**?

Delta-Xpress™ is a groundbreaking new splinting system that allows you to immobilise limbs in a fraction of the time needed, compared to conventional splints or casts.

The prefabricated design is intuitive to use, while providing secure and precise immobilisation.

Features and benefits

■ Precise clinical positioning

- Pre-defined angles hold limb in optimal healing position.
- Product comes ready-to-apply.
- Easy and secure application.

■ Quick – the fastest immobilisation ever

- Prefabricated design allows instant application.
- No activation, no curing, no clean-up.
- No need for water activation or even working space.
- Instructions and icons on product support quick and correct application.

■ Secure – safe and strong

- Rigid compound material.
- Good conformability for individual fit.
- Fixation with bandage or hook & loop-straps.

Indications where Delta-Xpress™ can be used:

- Post-traumatic emergency immobilisation.
- Immediate immobilisation before applying a long-term splint or cast.
- Primary care of soft tissue injuries.
- Post-operative immobilisation.



Splint designs

■ Wrist and finger splint

- Rigid stay for optimal wrist, hand and finger position.
- The stay has an angle of 100° to support the flexion of the MP joint and stabilises the wrist in 30° dorsal extension.



■ Wrist splint

- Rigid stay for optimal wrist position.
- Wrist angle in 30° dorsal extension guarantees a functional therapy position.



■ Lower leg splint

- Pre-adjusted 90° angle for ideal healing position of lower leg.
- Ankle and foot rest inlay for optimal foot position.
- Malleoli paddings for increased comfort.



■ Elbow splint

- Pre-adjusted 90° angle assures optimal healing position of upper arm, elbow and lower arm.



High performance design features

In addition to ultra quick, precise and secure immobilisation, this unique, patented product provides various additional benefits:

- **Cushioned** - padding inside for comfort and absorption
- **Breathability** - to help avoid skin maceration
- **Size and shape** - adjustable with standard scissors
- **Convenient** - splints come in two or three sizes with measuring scale and colour code on product
- **Easy access to wounds**