

# **Gypsona® S**

## **Plaster of Paris Bandage**



- Fast immobilisation of fractures
- High mass weight resulting in low plaster loss
- Strong and durable casts

# Gypsona® S

## Plaster of Paris Bandage

- Fast immobilisation of fractures
- High mass weight resulting in low plaster loss
- Strong and durable casts

GYPSONA® S Plaster of Paris Bandage is a fast setting material. It consists of a cotton gauze which is spread with gypsum and spooled on a round core. The bandages wet rapidly and evenly with no dry spots. They are easy to apply using established GYPSONA® bandage techniques. The smooth creamy plaster consistency produces a cast of even surface and excellent finish.

### Advantages

The quick setting time of GYPSONA® S bandage ensures fast immobilisation of fractures. This makes it ideal for use in busy fracture clinics and plaster rooms. The high mass weight of GYPSONA® S Bandage ensures a strong and durable cast, resulting in low plaster loss when the bandage is immersed.

### Indications

1. Splinting of all types of fractures: acute and secondary conditions
2. Immobilisation of diseased bones and joints
3. Corrections of deformities
4. Negative and positive casting techniques
5. Emergency splinting

### Presentation

Two GYPSONA® S bandages are packed in a clear polythene film, which protects them from humidity and water, and are then packed in a durable carton. This maintains the product in the best conditions until its use.



### AVAILABILITY

Size	Unit of sale
5cm x 2m	80 bandages
7.5cm x 2m	80 bandages
10cm x 2m	40 bandages
12.5cm x 2m	40 bandages
15cm x 2m	40 bandages
20cm x 2m	20 bandages
5cm x 2.7m	80 bandages
7.5cm x 2.7m	80 bandages
10cm x 2.7m	40 bandages
12.5cm x 2.7m	40 bandages
15cm x 2.7m	40 bandages
20cm x 2.7m	20 bandages
10cm x 20m	1 slab dispenser
12.5cm x 20m	1 slab dispenser
15cm x 20m	1 slab dispenser
20cm x 20m	1 slab dispenser

To receive further information, please complete the Inquiry Form designed to help us meet your requirements.