

# Cutimed<sup>®</sup> Cavity



The wound dressing  
for moderate to  
high exudate levels.



- atraumatic
- non-adhesive
- soft and conformable
- easy to apply
- highly absorbent
- easy to remove

# Cutimed® Cavity

## For moderate to high exudate levels.

Cutimed® Cavity is a sterile dressing, made from hydrophilic, absorbent polyurethane foam. It is soft, highly absorbent and conformable, creating the ideal environment for moist wound healing.



For use as wound cover.



For use as a cavity dressing.


### Indications

For the management of wounds including cavity wounds healing by secondary intention with moderate to high exudate including:

- Leg ulcers
- Diabetic ulcers
- Pressure ulcers
- Cavity wounds
- Excisions and dehisced surgical wounds

Cutimed® Cavity may be used under medical attention in the management of infected wounds.

### Tensile strength after fluid uptake



- One piece, atraumatic removal of Cutimed® Cavity

Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.

### Ordering information

REF No.	Size	Dressings/Box	HCPCS
7262100	2 x 2 3/8 in.	10	A6215
7262101	4 x 4 in.	10	A6215
7262103	6 x 6 in.	10	A6215

HCPCS Codes listed are only suggested codes; other codes may be more applicable depending on intended product use. Codes are subject to change without notice.

Please refer to the product label and / or package insert for full instructions on the safe use of these products.