

JOBST® LymphCARE



New!

Bandage set for
oedema management

Product information

JOBST® LymphCARE

Bandage set for oedema management



Why compression therapy with JOBST® LymphCARE?

Oedema of lymphatic or venous origin can be mobilised by external pressure. JOBST® LymphCARE is a multiple layer bandage system exerting low resting but high working pressure.

JOBST® LymphCARE contains all the necessary materials (including different compression bandages, padding and fixation material) that are needed for the professional bandaging of a limb (arm or leg) according to the principle of Complex Decongestive Physiotherapy (CDP). The single components of JOBST® LymphCARE have been reliably used in oedema management.



Indication

JOBST® LymphCARE is indicated for the management of lymphoedema or oedema of upper and/or lower extremities.

Product	Width	No. of components in JOBST® LymphCARE		
		Arm (Ref. No. 78639)	Leg (Ref. No. 78642)	Lower Leg (Ref. No. 78641)
Tricofix® D5 Elastic tubular bandage	6 cm	1	-	-
	8 cm	-	-	1
	10 cm	-	1	-
Elastomull® Elastic fixation bandage	4 cm	20	20	20
Artiflex® Synthetic padding bandage	10 cm	2	2	3
	15 cm	2	2	-
Comprilan® Textile elastic short stretch bandage	6 cm	1	1	1
	8 cm	1	1	1
	10 cm	5	3	3
	12 cm	-	4	-
Leukoplast® hospital Surgical tape	1.25 cm	1	1	1
JOBST® Foam Rubber pad	1 cm	1	1	1
Cloth bag	42 cm	1	1	1
Instructions for Use	-	1	1	1

A complete solution for compression bandaging of arm, leg and lower leg

Meeting the needs completely

All materials necessary for compression bandaging of either arm, leg or lower leg are compiled in JOBST® LymphCARE.

Distinctly arranged and ready for use

All products are arranged in a clear manner for ease of use.

Reliable quality

JOBST® LymphCARE has been developed in co-operation with a well-known lymphoedema clinic.

Components

Each JOBST® LymphCARE set contains a specific amount of the following products depending on the limb bandaged (arm, leg or lower leg):

Tricofix®

Elastic tubular bandage



Elastomull®

Elastic fixation bandage



Artiflex®

Synthetic padding bandage



Comprilan®

Textile elastic short stretch bandage



Leukoplast® hospital

Surgical tape



JOBST® foam

Rubber pad



Contraindications

Ischaemia, uncontrolled congestive heart failure, untreated septic phlebitis, phlegmasia coerulea dolens.

Cautions

Immobility (confinement to bed), skin infections, weeping dermatoses, incompatibility to fabric, impaired sensitivity of the limb. Leukoplast® hospital contains dry natural rubber.

Other recommendations

Incorrectly applied bandages may cause problems during wearing. Open wounds must be covered with an appropriate dressing before the bandages are put on.